


CATALOGO DETECTOR FLAT PANEL CARERAY 1500 Cwe (35 CM X 43 CM)



1500Cwe

Wireless Cassette

14" x 17" Imaging Area

Scintillator Matrix: 3072 x 2560

Pixel Pitch: 140 μm

DQE Performance: > 71% @ 0 lp/mm

Product Dimensions: 46.0 x 38.4 x 1.5 cm

Applications

- Medical and Veterinary
- Mobile Healthcare

Special Features

- Magnetic folding handle for easier transport
- Full-field AED for broad system compatibility
- Pre-calibrated panels for easy out-of-box imaging