


CATALOGO DETECTOR FLAT PANEL CARERAY 1800 Cwe (43 CM X 43 CM)



1800 Cwe

Wireless Cassette

17" x 17" Imaging Area

Scintillator Matrix: 3072 x 3072

Pixel Pitch: 140 μm

DQE Performance: > 70% @ 0 lp/mm

Product Dimensions: 46.0 x 46.0 x 1.5 cm

Applications

- Medical and Veterinary
- Mobile Healthcare

Special Features

- Magnetic folding handle for easier transport
- Full-field AED for broad system compatibility
- Dual band wireless – 2.4 and 5 GHz channels