


CATALOGO DETECTOR FLAT PANEL CARERAY 750 CW (24 CM X 30 CM)



750CW

Wireless Cassette

10" x 12" Imaging Area

Scintillator Matrix: 2048 x 2560

Pixel Pitch: 120 μm

DQE Performance: > 70% @ 0 lp/mm

Product Dimensions: 33.3 x 38.2 x 1.5 cm

Applications

- Medical – Exceptionally low dose for neonatal care
- Veterinary – Popular for field equine imaging

Special Features

- Full-Field AED for broad system compatibility
- Direct-Deposit CsI for easier imaging and broad system compatibility
- Built to IP64 dustproof and water-resistant standards