

# CATALOGO DETECTOR FLAT PANEL MAMMOGRAFICO CARERAY 750MC (24 CM X 30 CM)



750Mc

---

## Mammography

10" x 12" Imaging Area

Scintillator Matrix: 3072 x 4096

Pixel Pitch: 77  $\mu\text{m}$

## Applications

- Mammography
- Retrofit

## Special Features

- Slender cassette fits into standard bucky trays for easy system upgrades
- Supports Automatic Exposure Control for easier, more reliable results
- High Dynamic Range and direct-deposit CsI for accurate low dose imaging