

Portable Generators


ULTRA 100

- 5.0kW Powerful X-ray output (110kV / 100mA)
- **Auto SID measuring by Ultrasonic sensor**
- Easy access to DR solution
- Angle measuring by goniometer
- Reversible LED display
- Durable Hand-switch made of metal connector



ULTRA 100 HF

- 4.0kW Powerful X-ray output (110kV / 100mA)
- Light weight
- Reversible LED Display
- Easy access to DR Solution
- 110 - 240VAC Universal voltage





Trust Begins with
ECOVIEV

Technical Specifications

Category	Specifications	
	ULTRA 100	ULTRA 100HF
Output Power	5kW	4kW
Line Voltage	200~240VAC 50/60Hz, single phase	110 – 240VAC 50/60Hz single phase
Focal spot	1.8 mm	1.4 mm
Anode heat storage	28KJ	33KJ
kV range	40 ~ 110kV	40 ~ 110kV
mA range	25 ~ 100mA	20 ~ 100mA
mAs range	0.32 ~ 200mAs	0.32 ~ 100mAs
Collimator Type	Manual dual slit type	
Filtration	2.8 mm Al equivalent	2.77 mm Al equivalent
Size (mm)	396 x 295 x 204(H) mm	375 x 247 x 208(H) mm
Net weight	18 kg	14 kg

* Specifications are subject to change.



Standard Accessory



Aluminum Carrying Case